



Product Information

CAS Registry No. 70572-93-9

Product Name: **Ranatensin R**

Catalog No. **R0251**

Chemical Name:

Synonym(s):

Mol. formula: $C_{90}H_{134}N_{30}O_{24}S$

Formula wt: 2052.3

M.p.

B.p.

Purity: $\geq 98\%$

Solubility:

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

A peptide isolated from the skin of frogs.

Ser-Asn-Thr-Ala-Leu-Arg-Arg-
Tyr-Asn-Gln-Trp-Ala-Thr-Gly-
His-Phe-Met-NH₂

References:

GenScript Corporation. <http://www.genscript.com> (accessed June, 2006).

PRICING

1 mg	\$144.00
2 mg	\$244.80
5 mg	\$432.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.