



Product Information

CAS Registry No. 501-36-0

Product Name: Resveratrol

Catalog No. R1776

Chemical Name: 1,3-Benzenediol, 5-(2-(4-hydroxyphenyl)ethenyl)-, (E)-

Synonym(s): 3,4',5-Stilbenetrio;
3,4',5-Trihydroxystilbene

Mol. formula: C₁₄H₁₂O₃

Formula wt: 228.24

M.p. 254°C


B.p.

Purity: ≥98%

Solubility: Soluble in DMSO and ethanol.

Stability (in original shipping container):

Storage temp: -20 deg C



Product Description:

An antioxidant found in grapes, shown to have significant cancer chemopreventive and anti-tumor activity.

References:

Jang, M., Cai, L., Udeani, G.O., *et al.* *Science*. **275**:218-220 (1997).
Carbo, N., Costelli, P., Baccino, F.M, Lopez-Soriano, F.J.
Biochem. Ciophys. Res. Commun. **254**:739-43 (1999).
Gautam, S.C., Xu, Y.X., Dumaguin, M., *et al.* *Bone Marrow Transplant.* **25**:639-45 (2000).

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.

PRICING

100 mg	\$57.40
500 mg	\$241.20

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com