



Product Information

CAS Registry No. 4759-48-2

Product Name: 13-cis-Retinoic acid

Catalog No. R1779

Chemical Name: 3,7-Dimethyl-9-(2,6,6-trimethyl-1-cyclohexen-1-yl)-2-cis-4-trans-6-trans-8-trans-nonatetraenoic acid

Synonym(s): Isotretinoin; Isotrex; 13-cis-Vitamin A acid; Teriosal

Mol. formula: C₂₀H₂₈O₂

Formula wt: 300.44

M.p. 174°C-175°C


B.p.

Purity: ≥98%

Solubility: Sparingly soluble in water.

Stability (in original shipping container):

Storage temp: -20 deg C



Product Description:

A vitamin A analog that inhibits cell proliferation and induces cell differentiation.

References:

Kelloff, G.J., Crowell, J.A., Hawk, E.T., *et al.* *J. Cell. Biochem. Suppl.* **26**:54-71 (1996).

PRICING

100 mg	\$47.50
250 mg	\$90.30
500 mg	\$158.10

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.