



Product Information

CAS Registry No. 90098-04-7

Product Name: **Rebamipide**

Catalog No. **R1806**

Chemical Name: α -[(4-Chlorobenzoyl)-amino]-1,2-dihydro-2-oxo-4-quinolinepropanoic acid

Synonym(s): Proamipide; OPC-12759; Mucosta

Mol. formula: $C_{19}H_{15}ClN_2O_4$

Formula wt: 370.79

M.p. 288°C-290°C(dec)

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in dimethyl- formamide. Slightly soluble in methanol and ethanol.

Stability (in original shipping container):

Storage temp:

Product Description:

An antiulcer agent that enhances mucosal resistance. It increases endogenous prostaglandin gastric mucosa, and also scavenges oxygen-derived free radicals and inhibits their production.



References:

Otsubo K, Morita S, Uchida M et al. Chem Pharm Bull. 39:2906-9 (1991).

Iinuma S, Naito Y, Yoshikawa T et al. Dig Dis Sci. 43:35S-39S (1998).

PRICING

| | |
|------|----------|
| 1 g | \$30.80 |
| 5 g | \$123.20 |
| 25 g | \$431.20 |

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.