



Product Information

CAS Registry No. 11103-57-4

Product Name: all-trans-Retinol / Vitamin A

Catalog No. R1876

Chemical Name: 3,7-Dimethyl-9-(2,6,6-trimethyl-1-cyclohexen-1-yl)-2,4,6,8-nonatetraen-1-ol

Synonym(s): Alphasterol; Retinol; Avitol

Mol. formula: C₂₀H₃₀O

Formula wt: 286.45

M.p. 63°C-64°C

B.p.

Purity: ≥98%


Solubility: Soluble in alcohol, methanol, chloroform, ether, fats or oils.

Stability (in original shipping container):

Storage temp: +4 deg C

Product Description:

It inhibits cell proliferation and induces cell differentiation.



References:

Kelloff, G.J., Boone, C.W., Crowell, J.A., *et al.* *Cancer Epidemiol. Biomarkers Prev.* 3:85-98 (1994).

PRICING

100 mg	\$50.00
250 mg	\$90.00
500 mg	\$130.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.