



# Product Information

**CAS Registry No.** 127-47-9

**Product Name:** Retinyl acetate

**Catalog No.** R1878

**Chemical Name:**

**Synonym(s):** Vitamin A acetate; all-trans-Retinol acetate; Myvak

**Mol. formula:** C<sub>22</sub>H<sub>32</sub>O<sub>2</sub>

**Formula wt:** 328.49

**M.p.** 57°C-58°C


**B.p.**

**Purity:** ≥98%

**Solubility:** Soluble in ethanol.

**Stability** (in original shipping container):

**Storage temp:** +4 deg C



## Product Description:

A vitamin A analog that inhibits cell proliferation and induces cell differentiation. A chemopreventive against rat mammary tumors.

## References:

Grubbs, C.J., Eto, I., Juliana, M.M., Hardin, J.M., Whitaker, L.M. *Anticancer Res.* **10**:661-666 (1990).

## PRICING

5 g	\$41.60
25 g	\$125.80
100 g	\$363.20

*Bulk quantity available upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.