



# Product Information

**CAS Registry No.**

**Product Name:** RFRP-1 (Human)

**Catalog No.** R2369

**Chemical Name:**

**Synonym(s):**

**Mol. formula:**  $C_{67}H_{101}N_{19}O_{14}S$

**Formula wt:** 1428.73

**M.p.**

**B.p.**

**Purity:**

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:**

**Product Description:**

H-Met-Pro-His-Ser-Phe-Ala-  
Asn-Leu-Pro-Leu-Arg-Phe-NH<sub>2</sub>

**References:**

## PRICING

0.5 mg	\$38.40
1 mg	\$65.28
2.5 mg	\$115.20

*Bulk quantity available  
upon request*

## **LKT Laboratories, Inc.**

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.