



Product Information

CAS Registry No. 93674-99-8

Product Name: RGDV

Catalog No. R2514

Chemical Name: Arginyl-glycyl-aspartyl-valine

Synonym(s): Rgdv tetrapeptide

Mol. formula: $C_{17}H_{31}N_7O_7$

Formula wt: 445.5

M.p.

B.p.

Purity: $\geq 98\%$


Solubility:

Stability (in original shipping container):

Storage temp: -20 deg C

Product Description:

The 4-residue sequence found at the cell-adhesive region of the peptide vitronectin.



References:

Lin HB, Garcia-Echeverria C, Asakura S et al. Biomaterials. 13:905-914 (1992).

PRICING

5 mg	\$32.00
10 mg	\$54.40
25 mg	\$96.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.