



# Product Information

**CAS Registry No.**

**Product Name:** R-G-E-S

**Catalog No.** R2516

**Chemical Name:**

**Synonym(s):**

**Mol. formula:**  $C_{16}H_{29}N_7O_8$

**Formula wt:** 447.45

**M.p.**

**B.p.**

**Purity:**

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:**

## Product Description:

A tetrapeptide structurally similar to RGDS. RGDS is a synthetic peptide that represents the fibroblast binding site of fibronectin. RGEs is often used in fibroblast binding experiments.

## References:

Pasula R, Wisniowski P, Martin WJ 2nd. Infect Immun. 70:1287-92 (2002).

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.

H-Arg-Gly-Glu-Ser-OH

## PRICING

0.5 mg	\$25.60
1 mg	\$43.52
2.5 mg	\$76.80

*Bulk quantity available upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com