



Product Information

CAS Registry No. 13292-46-1

Product Name: Rifampicin

Catalog No. R3220

Chemical Name: 3-[[[4-Methyl-1-piperazinyl]imino]methyl]rifamycin

Synonym(s): Rifampin; Abrifam; Eremfat; Rifa; Rifaldin; Rifoldin

Mol. formula: $C_{43}H_{58}N_4O_{12}$

Formula wt: 822.94

M.p. 183°C-188°C (dec)

B.p.

Purity: ≥98%


Solubility: Soluble in DMSO ethyl acetate, methanol and THF. Slightly soluble in water, acetone, and carbon

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

An antibacterial from 3-formyl rifamycin SV.



References:

Maggi, *et al.* *Chemotherapy*. **11**:285 (1966).

PRICING

1 g	\$26.30
5 g	\$87.90
25 g	\$329.30

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.