



# Product Information

**CAS Registry No.** 14897-39-3

**Product Name:** Rifamycin SV monosodium salt

**Catalog No.** R3222

**Chemical Name:**

**Synonym(s):** Rifamycin SV sodium salt;  
Rifamastene; Rifocin

**Mol. formula:**  $C_{37}H_{46}NNaO_{12}$

**Formula wt:** 719.75

**M.p.**

**B.p.**

**Purity:**  $\geq 98\%$


**Solubility:** Soluble in water (pH7.2).

**Stability** (in original shipping container):

**Storage temp:** Ambient

## Product Description:

Antimicrobial agent inhibits DNA polymerase.



## References:

Furesz S. Arioli V. Pallanza R. *Antimicrob Agents Chemother.* 5:770-7 (1965).

Frolova LY. Meldrays YA. Kochkina LL. *et al. Nucleic Acids Res.* 4:523-38 (1977).

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.

## PRICING

1 g	\$29.40
5 g	\$95.10

*Bulk quantity available upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com