



Product Information

CAS Registry No. 1501-84-4

Product Name: Rimantadine hydrochloride

Catalog No. R3249

Chemical Name: α -Methyltricyclo[3.3.1.1^{3,7}]decane-1-methanamine

Synonym(s): Flumadine; Meradan; Roflual

Mol. formula: C₁₂H₂₁N

Formula wt: 179.3

M.p. 373°C-375°C

B.p.

Purity: \geq 98%


Solubility: Soluble in water.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

Rimantadine hydrochloride is used for influenza treatment and prevention..



References:

Wills, *et al. Antimicrob. Ag. Chemother.* **31**:826.

Tsunoda, *et al. Antimicrob. Ag. Chemother.* **5**:553 (1965).

Keyser, *et al. Infect. Dis. Clin. N. Am.* **1**:459 (1987).

PRICING

25 mg \$60.80

50 mg \$108.70

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.