



Product Information

CAS Registry No. 54029-12-8

Product Name: Ricobendazole

Catalog No. R3310

Chemical Name: (5-propylsulfinyl)-1H-benzimidazol-2-yl)-methyl ester carbamic acid

Synonym(s): (+-)-Albendazole sulfoxide;
Albendazole oxide

Mol. formula: C₁₂H₁₅N₃O₃S

Formula wt: 281.33

M.p. 226°C-228°C (dec)

B.p.

Purity: ≥98%


Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

A low-soluble anthelmintic drug.



References:

Castillo JA, Palomo-Canales J, Garcia JJ et al. Drug Dev Ind Pharm. 25:1241-8 (1999).

PRICING

10 g	\$39.50
25 g	\$74.00
100 g	\$215.60

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.