



Product Information

CAS Registry No. 115436721

Product Name: Risedronate sodium

Catalog No. R3373

Chemical Name: [1-Hydroxy-2-(3-pyridinyl)ethylidene]bisphosphonic acid monosodium salt

Synonym(s): Risedronic acid monosodium salt;
Actonel; Optinate

Mol. formula: $C_7H_{10}NO_7P_2Na$

Formula wt: 305.10

M.p.

B.p.

Purity: $\geq 98\%$


Solubility: Soluble in water. Insoluble in organic solvents.

Stability (in original shipping container):

Storage temp: Ambient

Product Description:

A pyridinyl biphosphate which inhibits osteoclast-mediated bone resorption.



References:

Boiser S., Ferreras M., Peyruchaud O., Margetto S. *Cancer Res.* **60**:2949-2954 (2000).

Goa KL., Balfour JA. *Drugs Aging.* **13**:83-91 (1998).

PRICING

100 mg \$86.00

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.