



# Product Information

**CAS Registry No.** 61281-38-7

**Product Name:** R(+) Schisandrin A

**Catalog No.** S0830

**Chemical Name:** Dibenzo(a,c)cyclooctene, 5,6,7,8-tetrahydro-1,2,3,10,11,12-hexamethoxy-6,7-dimethyl-, stereoisomer

**Synonym(s):** R(+) Schizandrin A; R(+) deoxyschisandrin A; R(+) deoxyschizandrin A

**Mol. formula:** C<sub>24</sub>H<sub>32</sub>O<sub>6</sub>

**Formula wt:** 416.51

**M.p.** 117°C-118°C

**B.p.**

**Purity:** ≥98%


**Solubility:**

**Stability** (in original shipping container):

**Storage temp:** +4 deg C

**Product Description:**

[α]<sub>D</sub>: +110° (c=1.15, CHCl<sub>3</sub>)



**References:**

## PRICING

5 mg	\$71.50
10 mg	\$117.10

*Bulk quantity available  
upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.