




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	3211-76-5
Product ID	S1845
Product Name	L-(+)-Selenomethionine / 2-Amino-4-(methylseleno)butanoic acid
Chemical Name	Butanoic acid, 2-amino-4-(methylseleno)-, (S)-
Synonym	Selenium methionine; Selenomethionine
Formula	C ₅ H ₁₁ NO ₂ Se
Formula Wt.	196.11
Melting Point	265°C(dec.)
Boiling Point	
Purity	≥98%
Solubility	Soluble in water. Slightly
Stability	
Store Temp	-20 deg C
Ship Temp	
Product Description	An inhibitor of 2-acetylaminofluorene-induced hepatocarcinogenesis.
References	Mukherjee, B., Ghosh, S., Chatterjee, M. <i>J. Exp. Ther. Oncol.</i> 14 :209-217 (1996).



Pricing

10 mg	\$19.30
25 mg	\$39.60
100 mg	\$131.50

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.