




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	13190-97-1
Product ID	S5746
Product Name	Solanesol
Chemical Name	(2E,6E,10E,14E,18E,22E,26E,30E)-3,7,11,15,19,23,27,31,35-Nonamethyl-2,6,10,14,18,22,26,30,34-hexa-triacontanonaen-1-ol
Synonym	Nonaisoprenol
Formula	C ₄₅ H ₇₄ O
Formula Wt.	631.07
Melting Point	41.5-42.5°C
Boiling Point	
Purity	≥92%
Solubility	Soluble in organic solvents.
Stability	
Store Temp	-20 deg C
Ship Temp	
Product Description	A trisesquiterpene alcohol present in tobacco leaf. It is a possible precursor of polycyclic aromatic hydrocarbons.
References	Severson RF, Ellington JJ, Schlotzhauer PF et al. J Chromatogr. 139:269-82 (1977).



Pricing

50 mg	\$89.70
100 mg	\$138.00
500 mg	\$448.40

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.