



Product Information

CAS Registry No.

Product Name: Somatostatin-28 (1-12)

Catalog No. S5752

Chemical Name:

Synonym(s):

Mol. formula: $C_{49}H_{81}N_{17}O_{19}S_1$

Formula wt: 1244.36

M.p.

B.p.

Purity:

Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

H-Ser-Ala-Asn-Ser-Asn-Pro-
Ala-Met-Ala-Pro-Arg-Glu-OH

References:

PRICING

1 mg	\$108.80
2 mg	\$185.60
5 mg	\$326.40

*Bulk quantity available
upon request*

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.