




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	92950-40-8
Product ID	T0093
Product Name	2'-Acetyltaxol
Chemical Name	Benzenepropanoic acid, α -(acetyloxy)- β -(benzoyl-amino)-,6,12b-bis(acetyloxy)-12-(benzyloxy)- 2a,3,4,4a,5,6,9,10,11,12,12a,12b-dodecahydro-4,11- dihydroxy-4a,8,13,13-tetramethyl-5-oxo-7,11-
Synonym	Taxol 2'-acetate
Formula	C ₄₉ H ₅₃ NO ₁₅
Formula Wt.	895.95
Melting Point	
Boiling Point	
Purity	≥98%
Solubility	
Stability	
Store Temp	-20 deg C
Ship Temp	Ambient
Product Description	
References	



Pricing

5 mg	\$227.70
10 mg	\$353.20
25 mg	\$689.90

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.