




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	27548-93-2
Product ID	T0095
Product Name	Baccatin III
Chemical Name	7,11-Methano-5H-cyclodeca(3,4)benz(1,2-b)oxet-5-one-6,12b-bis(acetyloxy)-12-(benzyloxy)-1,2a,3,4,4a,6,- 9,10,11,12,12a,12b-dodecahydro-4,9,11-trihydroxy-4a,- 8,13,13-tetramethyl-,(2aR-
Synonym	
Formula	C ₃₁ H ₃₈ O ₁₁
Formula Wt.	586.63
Melting Point	229°C-234°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in DMSO or
Stability	
Store Temp	-20 deg C
Ship Temp	Ambient
Product Description	Baccatin III induces apoptosis.
References	Miller, M.C. 3rd., Johnson, K.R., Willingham, M.C., Fan, W. <i>Cancer Chemother. Pharmacol.</i> 44 :444-52 (1999).



Pricing

5 mg	\$162.70
10 mg	\$295.70
25 mg	\$591.40

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.