




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	106463-17-6
Product ID	T0251
Product Name	Tamsulosin Hydrochloride
Chemical Name	5-[(2R)-2-[[2-(2-Ethoxy-phenoxy)-ethyl]amino]propyl]-2-methoxybenzenesulfonamide hydrochloride
Synonym	Flomax; Harnal; Pradif
Formula	C ₂₀ H ₂₈ N ₂ O ₅ S.HCl
Formula Wt.	444.98
Melting Point	228°C-230°C
Boiling Point	
Purity	≥98%
Solubility	Slightly soluble in water.
Stability	
Store Temp	
Ship Temp	
Product Description	An α ₁ -adrenoceptor antagonist used to treat benign prostatic hypertrophy.
References	Takenaka T, Fujikura T, Honda K et al. Yakugaku Zasshi. 115:773-89 (1995).



Pricing

10 mg	\$44.80
25 mg	\$76.20
100 mg	\$196.00

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.