




**LKT Laboratories, Inc**  
2233 University Avenue West  
St. Paul, Minnesota 55114 USA  
Ph: 651.644.8424 Fax: 651.644.8357  
www.lktlabs.com info@lktlabs.com

# Product Information

<b>CAS No.</b>	89785-84-2
<b>Product ID</b>	T0299
<b>Product Name</b>	Tazobactam sodium
<b>Chemical Name</b>	(2 <i>S</i> ,3 <i>S</i> ,5 <i>R</i> )-3-Methyl-7-oxo-3-(1 <i>H</i> -1,2,3-triazol-1-ylmethyl)-4-thia-1-azabicyclo[3.2.0]heptane-2-carboxylic acid 4,4-dioxide sodium salt
<b>Synonym</b>	YTR-830; CL-307579
<b>Formula</b>	
<b>Formula Wt.</b>	322.28
<b>Melting Point</b>	>170°(dec.)
<b>Boiling Point</b>	
<b>Purity</b>	≥98% Sp. rotation:+135°~
<b>Solubility</b>	Soluble in water.
<b>Stability</b>	
<b>Store Temp</b>	Ambient
<b>Ship Temp</b>	Ambient
<b>Product Description</b>	A β-lactamase inhibitor with antibacterial properties.
<b>References</b>	Bryson HM. Brogden RN. <i>Drugs</i> . 47:506-35 (1994).



## Pricing

100 mg	\$57.50
500 mg	\$216.60
1 g	\$393.50

**Bulk quantities available upon request**

## LKT Laboratories, Inc.

Address: 2233 University Ave. W.  
St. Paul, MN 55114  
USA

Phone: 888-LKT-LABS (558-5227)  
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.