




**LKT Laboratories, Inc**  
2233 University Avenue West  
St. Paul, Minnesota 55114 USA  
Ph: 651.644.8424 Fax: 651.644.8357  
www.lktlabs.com info@lktlabs.com

# Product Information

<b>CAS No.</b>	70024-40-7
<b>Product ID</b>	T1670
<b>Product Name</b>	Terazosin Hydrochloride dihydrate
<b>Chemical Name</b>	1-(4-Amino-6,7-dimethoxy-2-quinazolinyl)-4-(tetrahydro-2-furoyl)piperazine monohydrochloride dihydrate
<b>Synonym</b>	Heitrin; Hytracin; Hytrin; Hytrinex; Itrin; Urodie; Vasocard
<b>Formula</b>	C <sub>19</sub> H <sub>25</sub> N <sub>5</sub> O <sub>4</sub> ·HCl·2H <sub>2</sub> O
<b>Formula Wt.</b>	459.93
<b>Melting Point</b>	271-274°C
<b>Boiling Point</b>	
<b>Purity</b>	
<b>Solubility</b>	Soluble in water.
<b>Stability</b>	
<b>Store Temp</b>	
<b>Ship Temp</b>	
<b>Product Description</b>	An α <sub>1</sub> -selective blocker. Its inhibitory effect on prostate tumor growth may be the result of antiangiogenic activity.
<b>References</b>	Kyncl JJ. J Clin Pharmacol. 33:878-83 (1993). Pan SL, Guh JH, Huang YW et al. J Urol. 169:724-9 (2003).



## Pricing

50 mg	\$69.00
250 mg	\$276.00
1 g	\$827.90

**Bulk quantities available upon request**

## LKT Laboratories, Inc.

Address: 2233 University Ave. W.  
St. Paul, MN 55114  
USA

Phone: 888-LKT-LABS (558-5227)  
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.