



# Product Information

**CAS Registry No.** 52232-67-4

**Product Name:** Teriparatide Acetate

**Catalog No.** T1675

**Chemical Name:**

**Synonym(s):**

**Mol. formula:**  $C_{181}H_{291}N_{55}O_{51}S_2$

**Formula wt:** 4117.72

**M.p.**

**B.p.**

**Purity:**

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:**

## Product Description:

A synthetic peptide that consists the 1-34 amino acid fragment of human parathyroid hormone. It stimulates new bone formation, repairs structural defects and reduces risks of fractures in postmenopausal women with severe osteoporosis.

Ser-Val-Ser-Glu-Ile-Gln-Leu-Met-His-Asn-Leu-Gly-Lys-His-Leu-Asn-Ser-Met-Glu-Arg-Val-Glu-Trp-Leu-Arg-Lys-Lys-Leu-Gln-Asp-Val-His-Asn-Phe-OH

## References:

Yodfat Y. Harefuah. 146:134-9, 164 (2007).

## PRICING

*Bulk quantity available upon request*

### LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.