




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	113712-98-4
Product ID	T1754
Product Name	Tenatoprazole
Chemical Name	5-Methoxy-2-[[[(4-Methoxy-3,5-Dimethyl-2-pyridinyl)methyl] sulphanyl]-1H-imidazo[4,5-b]pyridine
Synonym	Protop; Ulsacare
Formula	C ₁₆ H ₁₈ N ₄ O ₃ S
Formula Wt.	346.41
Melting Point	173°C-175°C
Boiling Point	
Purity	≥98%
Solubility	Soluble in dichlormethane.
Stability	
Store Temp	+4 deg C
Ship Temp	
Product Description	A proton pump inhibitor.
References	Nakamura T, Nippon Rinsho. 60 Suppl 2:650-4 (2002).



Pricing

100 mg	\$48.40
500 mg	\$144.90
1 g	\$231.80

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.