




LKT Laboratories, Inc
2233 University Avenue West
St. Paul, Minnesota 55114 USA
Ph: 651.644.8424 Fax: 651.644.8357
www.lktlabs.com info@lktlabs.com

Product Information

CAS No.	148-79-8
Product ID	T2930
Product Name	Thiabendazole
Chemical Name	2-(4-Thiazolyl)-1H-benzimidazole
Synonym	MK-360; Thibenzole; Equizole; Mertect; Storite; TBZ; Tecto.
Formula	C ₁₀ H ₇ N ₃ S
Formula Wt.	201.25
Melting Point	304°C-305 °C
Boiling Point	
Purity	≥98%
Solubility	Soluble in DMF and DMSO.
Stability	Stable at room temperature.
Store Temp	
Ship Temp	
Product Description	A widely used fungicide and antihelminthal agent. Mediates Cytochrome P450 1A1 induction.
References	Ogunsusi RA. Res Vet Sci. 25:251-252 (1978). Iemaire G, Delescluse C, Pralavorio M et al. Life Sci. 74:2265-2278 (2004).



Pricing

10 g	\$31.40
100 g	\$61.50
500 g	\$223.30

Bulk quantities available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W.
St. Paul, MN 55114
USA

Phone: 888-LKT-LABS (558-5227)
651-644-8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.