



# Product Information

**CAS Registry No.**

**Product Name:** Uperolein

**Catalog No.** U6118

**Chemical Name:**

**Synonym(s):**

**Mol. formula:**  $C_{57}H_{79}N_{13}O_{16}S$

**Formula wt:** 1234.42

**M.p.**

**B.p.**

**Purity:**

**Solubility:**

**Stability** (in original shipping container):

**Storage temp:**

**Product Description:**

A physalaemin-like endecapeptide that is selective for Neurokinin 1 receptor.

pGlu-Pro-Asp-Pro-Asn-Ala-  
Phe-Tyr-Gly-Leu-Met-NH<sub>2</sub>

**References:**

Dike A, Cowsik SM. J Struct Biol. 156: 442-52 (2006).

## PRICING

1 mg	\$70.40
2 mg	\$120.00
5 mg	\$211.20

*Bulk quantity available  
upon request*

## LKT Laboratories, Inc.

**Address:** 2233 University Ave. W  
St. Paul, MN 55114  
USA

**Phone:** 888-LKT-Labs (558-5227)  
651 644 8424

**Fax:** 651-644-8357

**Email:** info@lktlabs.com

**Web:** www.lktlabs.com

**Caution:** This product is intended for laboratory and research use only. It is not for human or drug use.