



Product Information

CAS Registry No.

Product Name: Urocortin, rat

Catalog No. U6855

Chemical Name:

Synonym(s):

Mol. formula: $C_{206}H_{338}N_{62}O_{64}$

Formula wt: 4707.37

M.p.

B.p.

Purity:

Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

H-Asp-Asp-Pro-Pro-Leu-Ser-Ile-Asp-Leu-Thr-Phe-His-Leu-Leu-Arg-Thr-Leu-Leu-Glu-Leu-Ala-Arg-Thr-Gln-Ser-Gln-Arg-Glu-Arg-Ala-Glu-Gln-Asn-Arg-Ile-Ile-Phe-Asp-Ser-Val-NH₂

References:

PRICING

0.5 mg	\$160.00
1 mg	\$272.00
2.5 mg	\$480.00

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.