



Product Information

CAS Registry No.

Product Name: Vesicular Stomatitis Virus peptide

Catalog No. V1872

Chemical Name:

Synonym(s):

Mol. formula: $C_{44}H_{66}N_{12}O_{12}$

Formula wt: 955.09

M.p.

B.p.

Purity:

Solubility:

Stability (in original shipping container):

Storage temp:

Product Description:

A antigenic peptide of vesicular stomatitis virus binds to H-2Kb antigens.

Arg-Gly-Tyr-Val-Tyr-Gln-Gly-Leu

References:

Shibata K, Imarai M, van Bleek GM et. al. Proc Natl Acad Sci U S A. 89: 3135-9 (1992).

PRICING

1 mg	\$32.00
2 mg	\$59.20
5 mg	\$105.60

Bulk quantity available upon request

LKT Laboratories, Inc.

Address: 2233 University Ave. W
St. Paul, MN 55114
USA

Phone: 888-LKT-Labs (558-5227)
651 644 8424

Fax: 651-644-8357

Email: info@lktlabs.com

Web: www.lktlabs.com

Caution: This product is intended for laboratory and research use only. It is not for human or drug use.